


## Technical Information about Rp-8-Br-cAMPS

Potent membrane-permeant, metabolically stable inhibitor of cAMP-dependent protein kinases

Update: June 08, 2017\_HU



### Abbreviation:

**Rp-8-Br-cAMPS**

| Formula   | CAS No.       | Molecular Weight | UV                                       | BIOLOG Cat.No. |
|---|---------------|------------------|--|----------------|
| C <sub>10</sub> H <sub>10</sub> BrN <sub>5</sub> O <sub>5</sub> PS·Na | [129735-00-8] | 446.2            | λ <sub>max</sub> 264 nm / ε 17000 / pH 7 | B 001          |

**Name:** 8- Bromoadenosine- 3', 5'- cyclic monophosphorothioate, Rp-isomer.

**Description:** Rp-8-Br-cAMPS is an analogue of the parent compound cyclic AMP in which the hydrogen in position 8 of the nucleobase is replaced by bromine and the equatorial one of the two exocyclic oxygen atoms in the cyclic phosphate moiety is modified by sulfur. The suffix "p" indicates that R/S nomenclature refers to phosphorus.

**Properties:** Rp-8-Br-cAMPS is a combination of the well accepted protein kinase inhibitor Rp-cAMPS with the widely used 8-bromo cyclic AMP, resulting in a membrane-permeant cyclic AMP antagonist which is not metabolized by mammalian cyclic nucleotide phosphodiesterases.

Rp-8-Br-cAMPS is about 1.5 times more lipophilic compared to 8-Br-cAMP, and 2 times and 3 times more compared to Rp-cAMPS and cAMP, respectively.

In contrast to common ATP-site inhibitors, Rp-8-Br-cAMPS discriminates between both isozymes of protein kinase A preferring type I (Gjertsen et al., 1995) and thus provides additional selectivity. By occupying cAMP binding sites Rp-8-Br-cAMPS prevents the kinase holoenzyme from dissociation and thus from activation.

Due to this working principle preincubation of the inhibitor prior to the activation step is necessary for optimal results.

**Application:** Experience shows that applicable concentrations of Rp-8-Br-cAMPS depend on the type of biosystem, its membrane properties and kinase content. A main application for Rp-8-Br-cAMPS is to eliminate the first messenger-stimulated phosphorylation by cyclic AMP-dependent protein kinase. For this purpose preincubation (e.g. 20 min.) is important, since the production of intracellular cyclic AMP initiated by a first messenger is much faster than the antagonist can penetrate the membrane when given extracellularly. Since Rp-8-Br-cAMPS is hydrolytically stable in mammalian and many other systems there is no danger of degradation during incubation periods. If you have good or moderate results with 8-Br-cAMP or Rp-cAMPS, you can be sure that Rp-8-Br-cAMPS will be membrane-permeant in your system as well.

**Specification:** Lyophilized or crystallized sodium salt. The free acid or other salt forms are available upon request. Equal concentrations of Rp-8-Br-cAMPS can appear very different in volume due to sensitivity of the lyophilized form to humidity. The compound can even contract to small volume droplets. Normally the product is located in the conical bottom of the tube. Micromolar quantities are determined by UV at λ<sub>max</sub>. BIOLOG LSI also offers the corresponding agonistic Sp- isomer (Sp-8-Br-cAMPS; Cat. No. B 002) and 8-Br-cAMP (Cat. No. B 007).


**Purity:** Typical analysis is better than 99% (HPLC / UV / 264 nm). The product has been specially treated for low endotoxin levels, however, we do not guarantee total absence of lipopolysaccharides. The content of the vial is not sterile.

**Caution:** Since even minor impurities of 8-Br-cAMP (0.2%) or Sp-8-Br-cAMPS can already activate protein kinase A and compete with the antagonistic effect of the Rp-isomer, it is very important to work with strictly pure compounds concerning cyclic nucleotide contaminants. Rp-8-Br-cAMPS is specially checked for absence of both, Sp-8-Br-cAMPS and 8-Br-cAMP (< 0.05% when packed). However, we cannot guarantee total absence of 8-Br-cAMP due to its formation during prolonged storage.

**Stability and Storage:** Rp-8-Br-cAMPS has sufficient stability at room temperature and does not need special care during handling or shipment. Nevertheless, we recommend that the compound should be stored in the freezer, for longer storage periods preferably in freeze-dried form, since desulfurization yielding 8-Br-cAMP can occur slowly. This normally rather slow sulfur/oxygen exchange is accelerated by oxidizing agents. Bright light, radioactivity or UV radiation should be avoided.

**Solubility:** Detailed information on the solubility of Rp-8-Br-cAMPS in water and various buffers are listed in the solubility chart below. Concentrations have been tested at ambient temperature and can be considered as minimum concentrations usually obtainable, however, slight batch-to-batch variations cannot be ruled out. When opening the tube please make sure that no substance is lost within the cap. Please rinse tube walls carefully and preferably use ultrasonic or vortex to achieve total and uniform mixing.

| No.  | Solvent  | Solubility [mM] |
|------|--|-----------------|
| I    | H <sub>2</sub> O                                 | 250             |
| II   | DMSO   | 250             |
| III  | DMF  | 250             |
| IV   | Ethanol 96%                                      | 250             |
| V    | Methanol   | 250             |
| VI   | PBS, pH 7.4                                      | 250             |
| VII  | 100 mM Na <sub>2</sub> HPO <sub>4</sub> , pH 7.0 | 250             |
| VIII | 25 mM Hepes/NaOH, pH 7.2                         | 250             |
| IX   | 25 mM Tris/HCl, pH 7.4                           | 250             |



**Toxicity and Safety:** Since cyclic AMP has multiple tasks in every organism it is very likely that lipophilic cAMP analogues will interfere with many cell regulation processes *in vivo*. However, due to the rather small quantities to work with no health hazards have been reported. Nevertheless, please keep in mind that the *in vivo* properties of these compounds are not sufficiently characterized up to now. Avoid skin contact or ingestion and allow only trained personnel to handle the product.

Our products are designed, developed and sold for research purposes only! They are intended for *in vitro* and nonhuman *in vivo* laboratory applications. Any other use requires approval of health authorities.

**Not for drug, household or related uses!**

#### Selected References for Rp-8-Br-cAMPS:

For a detailed list please inquire or visit our website (<http://www.biolog.de>).

Herfindal, L.; Nygaard, G.; Kopperud, R.; Krakstad, C.; Doeskeland, S.O.; Selheim, F., *Biochem. Biophys. Res. Commun.*, **437**, 603 - 608 (2013): "Off-target Effect of the Epac Agonist 8-pCPT-2'-O-Me-cAMP on P2Y<sub>12</sub> Receptors in Blood Platelets"

Chepurny, O.G.; Bertinetti, D.; Diskar, M.; Leech, C.A.; Afshari, P.; Tsalkova, T.; Cheng, X.; Schwede, F.; Genieser, H.-G.; Herberg, F.W.; Holz, G.G., *Mol. Endocrinol.*, **27**, 1267-1282 (2013): "Stimulation of Proglucagon Gene Expression by Human GPR119 Enteroendocrine L-Cell Line GLUTag"

Werner, K.; Schwede, F.; Genieser, H.-G.; Geiger, J.; Butt, E., *Naunyn Schmiedebergs Arch. Pharmacol.*, **384**, 169 - 176 (2011): "Quantification of cAMP and cGMP Analogs in Intact Cells: Pitfalls in Enzyme Immunoassays for Cyclic Nucleotides"

Rah, S.-Y.; Mushtaq, M.; Nam, T.-S.; Kim, S.H.; Kim, U.-H., *J. Biol. Chem.*, **285**, 21877 - 21887 (2010): "Generation of Cyclic ADP-ribose and Nicotinic Acid Adenine Dinucleotide Phosphate by CD38 for Ca<sup>2+</sup> Signaling in Interleukin-8-treated Lymphokine-activated Killer Cells"

Poppe, H.; Rybalkin, S.D.; Rehmann, H.; Hinds, T.R.; Tang, X.-B.; Christensen, A.E.; Schwede, F.; Genieser, H.-G.; Bos, J.L.; Doeskeland, S.O.; Beavo, J.A.; Butt, E., *Nature Methods*, **5**, 277 - 278 (2008): "Cyclic Nucleotide Analogs as Probes of Signaling Pathways"

Petersen, R.K.; Madsen, L.; Pedersen, L.M.; Hallenborg, P.; Hagland, H.; Viste, K.; Doeskeland, S.O.; Kristiansen, K., *Mol. Cell. Biol.*, **28**, 3804 - 3816 (2008): "Cyclic AMP (cAMP)-Mediated Stimulation of Adipocyte Differentiation Requires the Synergistic Action of Epac- and cAMP-Dependent Protein Kinase-Dependent Processes"

Ouyang, M.; Zhang, L.; Zhu, J.; Schwede, F.; Thomas, S.A., *PNAS*, **105**, 11993 - 11997 (2008): "Epac Signaling is Required for Hippocampus-dependent Memory Retrieval"

Bryn, T.; Mahic, M.; Aandahl, E.M.; Froland, S.S.; Aukrust, P.; Tasken, K., *AIDS Res. Hum. Retroviruses*, **24**, 1013 - 1015 (2008): "Inhibition of Protein Kinase A Improves Effector Function of Monocytes from HIV-Infected Patients"

Oestreich, E.A.; Wang, H.; Malik, S.; Kaproth-Joslin, K.A.; Blaxall, B.C.; Kelley, G.G.; Dirksen, R.T.; Smrcka, A.V., *Am. Soc. Biochem. Mol. Biol.* **2** (2006): **282** 5488-95 "Epac-mediated Activation of Phospholipase C-epsilon Plays a Critical Role in Beta-adrenergic Receptor Dependent Enhancement of Ca<sup>2+</sup> Mobilization in Cardiac Myocytes"

Kang, G.X.; Chepurny, O.G.; Malester, B.; Rindler, M.J.; Rehmann, H.; Bos, J.L.; Schwede, F.; Coetzee, W.A.; Holz, G.G., *J. Physiol.-London*, **573**, 595 - 609 (2006): "cAMP Sensor Epac as a Determinant of ATP-sensitive Potassium Channel Activity in Human Pancreatic Beta Cells and Rat INS-1 Cells"

Karege, F.; Schwald, M.; El Kouaissi, R., *Neuropsychopharmacol.*, **29**, 805 - 812 (2005): "Drug-induced Decrease of Protein Kinase A Activity in BDNF Expression of Bipolar Affective Disorder"

- Jensen, B.O.; Selheim, F.; Doeskeland, S.O.; Gear, A.R.L.; Holmsen, H., *Blood*, **104**, 2775-2782 (2004): "Protein Kinase A Mediates Inhibition of the Thrombin-Induced Platelet Shape Change by Nitric Oxide"
- Johansson, C.C.; Yndestad, A.; Enserink, J.M.; Ree, A.H.; Aukrust, P.; Tasken, K., *Endocrinology*, **145**, 5177 - 5184 (2004): "The Epidermal Growth Factor-like Growth Factor Amphiregulin is Strongly Induced by the Adenosine 3',5'-monophosphate Pathway in Various Cell Types"
- Horrigan, L.A.; Kelly, J.P.; Connor, T.J., *Int. Immunopharmacology*, **4**, 1409 - 1417 (2004): "Caffeine Suppresses TNF-alpha Production via Activation of the cyclic AMP/Protein Kinase A Pathway"
- Rodríguez-Moreno, A.; Sihra, T.S., *J. Physiol.-London*, **557**, 733 - 745 (2004): "Presynaptic Kainate Receptor Facilitation of Glutamate Release Involves Protein Kinase A in the Rat Hippocampus"
- Lu, Y.; Li, Y.; Herin, G.A.; Aizenman, E.; Epstein, P.M.; Rosenberg, P.A., *Eur. J. Neurosci.*, **19**, 2669 - 2681 (2004): "Elevation of Intracellular cAMP Evokes Activity-dependent Release of Adenosine in Cultured Rat Forebrain Neurons"
- Jung, S.R.; Kim, M.H.; Hille, B.; Nguyen, T.D.; Koh, D.S., *Am. J. Physiol.-Cell Physiology*, **286**, C573 - C579 (2004): "Regulation of Exocytosis by Purinergic Receptors in Pancreatic Epithelial Cells"
- Tiwari, S.; Felekakis, K.; Moon, E.Y.; Flies, A.; Sherr, D.H.; Lerner, A., *Blood*, **103**, 2661 - 2667 (2004): "Among Circulating Hematopoietic Cells, B-CLL Uniquely Expresses Functional EPAC1 but EPAC1-mediated Rap1 Activation does not Account for PDE4 Inhibitor-induced Apoptosis"
- Kageyama, K.; Furukawa, K.I.; Miki, I.; Terui, K.; Motomura, S.; Suda, T., *J. Cardiovasc. Pharmacol.*, **42**, 561 - 565 (2003): "Vasodilative Effects of Urocortin II via Protein Kinase A and a Mitogen-activated Protein Kinase in Rat Thoracic Aorta"
- Eliasson, L.; Ma, X.S.; Renström, E.; Barg, S.; Berggren, P.O.; Galvanovskis J.; Gromada, J.; Jing, X.J.; Lundquist, I.; Salehi, A.; Sewing, S.; Rorsman P., *J. Gen. Physiol.*, **121**, 181 - 197 (2003): "SUR1 Regulates PKA-Independent cAMP-Induced Granule Priming in Mouse Pancreatic B-Cells"
- Pacheco, M.A.; Pastoor, T.E.; Wecker, L., *Mol. Brain. Res.*, **114**, 65 - 72 (2003): "Phosphorylation of the Alpha 4 Subunit of Human Alpha 4 Beta(2) Nicotinic Receptors: Role of cAMP-dependent Protein Kinase (PKA) and Protein Kinase C (PKC)"
- Schmidtko, A.; Ruth, P.; Geisslinger, G.; Tegeder, I., *Nitric Oxide-Biol. Chem.*, **8**, 89 - 94 (2003): "Inhibition of Cyclic Guanosine 5'-Monophosphate-dependent Protein Kinase I (PKG-I) in Lumbar Spinal Cord Reduces Formalin-induced Hyperalgesia and PKG Upregulation"
- Rozenfeld, C.; Martinez, R.; Figueiredo, R.T.; Bozza, M.T.; Lima, F.R.S.; Pires, A.L.; Silva, P.M.; Bonomo, A.; Lannes-Vieira, J.; De Souza, W.; Moura-Neto, V., *Infect. Immun.*, **71**, 2047 - 2057 (2003): "Soluble Factors Released by *Toxoplasma gondii*-infected Astrocytes Down-modulate Nitric Oxide Production by Gamma Interferon-activated Microglia and Prevent Neuronal Degeneration"
- Lin, S.L.; Johnson-Farley, N.N.; Lubinsky, D.R.; Cowen, D.S., *J. Neurochem.*, **87**, 1076 - 1085 (2003): "Coupling of Neuronal 5-HT<sub>7</sub> Receptors to Activation of Extracellular-regulated Kinase through a Protein Kinase A-independent Pathway that can Utilize Epac"
- Lu, Z.; Kolodecik, T.R.; Karne, S.; Nyce, M.; Gorelick, F., *Am. J. Physiol. Gastrointest. Liver Physiol.*, **285**, G822 - G828 (2003): "Effect of Ligands that Increase cAMP on Caerulein-induced Zymogen Activation in Pancreatic Acini"
- Hareng, L.; Meergans, T.; von Aulock, S.; Volk, H. D.; Hartung, T., *Eur. J. Immunol.*, **33**, 2287 - 2296 (2003): "Cyclic AMP Increases Endogenous Granulocyte Colony-Stimulating Factor Formation in Monocytes and THP-1 Macrophages Despite Attenuated TNF-Alpha Formation"
- Holm, A.M.; Aukrust, P.; Aandahl, E.M.; Müller, F.; Tasken, K.; Frøland, S.S., *J. Immunol.*, **170**, 5772 - 5777 (2003): "Impaired Secretion of IL-10 by T Cells from Patients with Common Variable Immunodeficiency-Involvement of Protein Kinase A Type I"
- Gendron, L.; Oligny, J.F.; Payet, M.D.; Gallo-Payet, N., *J. Biol. Chem.*, **278**, 3606 - 3614 (2003): "Cyclic AMP-independent Involvement of Rap1/B-Raf in the Angiotensin II AT(2) Receptor Signaling Pathway in NG108-15 Cells"
- Chepurny, O.G.; Hussain, M.A.; Holz, G.G., *Endocrinology*, **143**, 2303 - 2313 (2002): "Exendin-4 as a Stimulator of Rat Insulin I Gene Promoter Activity via bZIP/CRE Interactions Sensitive to Serine/Threonine Protein Kinase Inhibitor Ro 31-8220"
- Sasaki, S.; Nakagaki, I.; Kondo, H.; Hori, S., *Pflugers Arch.*, **445**, 342 - 351 (2002): "Involvement of the Ryanodine-sensitive Ca<sup>2+</sup> Store in GLP-1-induced Ca<sup>2+</sup> Oscillations in Insulin-secreting HIT Cells"
- Matuszyk, J.; Cebrat, M.; Kalas, W.; Strzadala, L., *Int. Immunopharmacol.*, **2**, 435 - 442 (2002): "HA1004, an Inhibitor of Serine/Threonine Protein Kinases, Restores the Sensitivity of Thymic Lymphomas to Ca<sup>2+</sup>-mediated Apoptosis through a Protein Kinase A-independent Mechanism"
- Bagley, K.C.; Abdelwahab, S.F.; Tuskan, R.G.; Fouts, T.R.; Lewis, G.K., *Infect. Immun.*, **70**, 5533 - 5539 (2002): "Cholera Toxin and Heat-Labile Enterotoxin Activate Human Monocyte-derived Dendritic Cells and Dominantly Inhibit Cytokine Production through a Cyclic AMP-dependent Pathway"
- Christensen, A.E.; de Rooij, J.; Kopperud, R.; Nilsen, O.K.; Dao, K.K.; Schwede, F.; Bos, J.L.; Genieser, H.-G.; Døskeland, S.O., Abstracts of the 13th Protein Kinase Symposium, Oslo/Norway, 2002: "Comparison of the cAMP Binding Sites of EPAC-1 and cAMP-Kinase"
- Johansson, C.C.; Bryn, T.; Yndestad, A.; Eiken, H.G.; Bjerkeli, V.; Frøland, S.S.; Aukrust, P.; Tasken, K., Abstracts of the 13th Protein Kinase Symposium, Oslo/Norway, 2002: "Cytokine Networks are Pre-activated in T Cells from HIV-infected Patients and under Control of cAMP as Assessed by Cytokine Arrays"
- Johnston, C.A.; Beazely, M.A.; Vancura, A.F.; Wang, J.K.T.; Watts, V.J., *J. Neurochem.*, **82**, 1087 - 1096 (2002): "Heterologous Sensitization of Adenylate Cyclase is Protein Kinase A-dependent in Cath.a Differentiated (CAD)-D-2L Cells"

- Aandahl, E.M.; Moretto, W.J.; Haslett, P.A.; Vang, T.; Bryn, T.; Tasken, K.; Nixon, D.F., *J. Immunol.*, **169**, 802 - 808 (2002): "Inhibition of Antigen-specific T Cell Proliferation and Cytokine Production by Protein Kinase A Type I"
- Martin, M.C.; Dransfield, I.; Haslett, C.; Rossi, A.G., *J. Biol. Chem.*, **276**, 45041 - 45050 (2001): "Cyclic AMP Regulation of Neutrophil Apoptosis Occurs via a Novel Protein Kinase A-independent Signaling Pathway"
- Braun, K.W.; Tribley, W.A.; Griswold, M.D.; Kim, K.H., *J. Biol. Chem.*, **275**, 4145 - 4151 (2000): "Follicle-stimulating Hormone Inhibits All-trans-retinoic Acid-induced Retinoic Acid Receptor  $\alpha$  Nuclear Localization and Transcriptional Activation in Mouse Sertoli Cell Lines"
- Spicuzza, L.; Belvisi, M.G.; Birrell, M.A.; Barnes, P.J.; Hele, D.J.; Giembycz, M.A., *Br. J. Pharmacol.*, **133**, 1201 - 1212 (2001): "Evidence that the Anti-Spasmogenic Effect of the Beta-Adrenoceptor Agonist, Isoprenaline, on Guinea-Pig Trachealis is not Mediated by Cyclic AMP-dependent Protein Kinase"
- Zheng, J.; Ghorpade, A.; Niemann, D.; Cotter, R.L.; Thylin, M.R.; Epstein, L.; Swartz, J.M.; Shepard, R.B.; Liu, X.; Nukuna, A.; Howard, E., *J. Virol.*, **73**, 8256 - 8267 (1999): "Lymphotropic Virions Affect Chemokine Receptor-Mediated Neural Signaling and Apoptosis: Implications for Human Immunodeficiency Virus Type 1-Associated Dementia"
- Kurosaki, K.; Ikeda, U.; Maeda, Y.; Shimpo, M.; Ueno, S.; Shimada, K., *Cardiovasc. Res.*, **38**, 192 - 197 (1998): "Effects of Vesnarinone on Nitric Oxide Synthesis in Rat Cardiac Myocytes"
- Keryer, G.; Alsat, E.; Tasken, K.; Evain-Brion, D., *J. Cell Sci.*, **111**, 995 - 1004 (1998): "Cyclic AMP-dependent Protein Kinases and Human Trophoblast Cell Differentiation in Vitro"
- Aandahl, E. M.; Aukrust, Pal, Skalhegg, B.S.; Müller, F.; Froland, S.S.; Hansson, V.; Tasken, K., *FASEB J.*, **12**, 855 - 862 (1998): "Protein Kinase A Type I Antagonists Restores Immune Responses of T Cells from HIV-infected Patients"
- Brar, A.K.; Frank, G.R.; Kessler, C.A.; Cedars, M.I.; Handwerger, S., *Endocrine* **6**, 301 - 307 (1997): "Progesterone-dependent Decidualization of the Human Endometrium is Mediated by cAMP"
- Derkinderen, P.; Toutant, M.; Burgaya, F.; Le Bert, M.; Siciliano, J.C.; de Franciscis, V.; Gelman, M.; Girault, J.-A., *Science*, **273**, 1719 - 1722 (1996): "Regulation of a Neuronal Form of Focal Adhesion Kinase by Anandamide"
- Hemmings, H.C. in *Regulatory Protein Modification: Techniques and Protocols; Neuromethods* **30**, 121 - 218 (1996); H.C. Hemmings (ed.), Humana Press Inc.: "Protein Kinase and Phosphatase Inhibitors. Applications in Neuroscience"
- Wilderman, M.J.; Armstead, W.M., *Eur. J. Pharmacol.*, **309**, 243 - 249 (1996): "Influence of cAMP on Cerebrospinal Fluid Opioid Concentration: Role in cAMP-Induced Pial Artery Dilation"
- Stella, N.; Magistretti, P.J., *J. Biol. Chem.*, **271**, 23705 - 23710 (1996): "Vasoactive Intestinal Peptide (VIP) and Pituitary Adenylate Cyclase-activating Polypeptide (PACAP) Potentiate the Glutamate-evoked Release of Arachidonic Acid from Mouse Cortical Neurons- Evidence for a cAMP-independent Mechanism"
- Murphy-Ullrich, J. E.; Pallero, M.A.; Boerth, N.; Greenwood, J.A.; Lincoln, T.M.; Cornwell, T.L., *J. Cell Sci.*, **109**, 2499 - 2508 (1996): "Cyclic GMP-dependent Protein Kinase is Required for Thrombospondin and Tenascin Mediated Focal Adhesion Disassembly"
- Rooney, T.A.; Joseph, S.K.; Queen, C.; Thomas, A.P., *J. Biol. Chem.*, **271**, 19817 - 19825 (1996): "Cyclic GMP Induces Oscillatory Calcium Signals in Rat Hepatocytes"
- Ikeda, U.; Ikeda, M.; Kano, S.; Kanbe, T.; Shimada, K., *Eur. J. Pharmacol.*, **314**, 197 - 202 (1996): "Effect of Cilostazol, A cAMP Phosphodiesterase Inhibitor, on Nitric Oxide Production by Vascular Smooth Muscle Cells"
- Holen, I.; Gordon, P.B.; Stroemhaug, P.E.; Seglen, P.O., *Eur. J. Biochem.*, **236**, 163 - 170 (1996): "Role of cAMP in the Regulation of Hepatocytic Autophagy"
- Armstead, W.M., *Am. J. Physiol.*, **40**, H166 - H172 (1996): "cGMP and cAMP in Prostaglandin-induced Pial Artery Dilation and Increased CSF Opioid Concentration"
- Wolf, F.I.; diFrancesco, A.; Covacci, V.; Corda, D.; Cittadini, A., *Arch. Biochem. Biophys.*, **331**, 194 - 200 (1996): "Regulation of Intracellular Magnesium in Ascites Cells: Involvement of Different Regulatory Pathways"
- Tenor, H.; Hatzelmann, A.; Church, M.K.; Schudt, C.; Shute, J.K., *Br. J. Pharmacol.*, **118**, 1727 - 1735 (1996): "Effects of Theophylline and Rolipram on Leukotriene C4 (LTC4) Synthesis and Chemotaxis of Human Eosinophils from Normal and Atopic Subjects"
- Cho-Chung, Y.S.; Jastorff, B.; Genieser, H.-G., USA patent application, 1992: "Phosphorothioate Derivatives of Cyclic AMP Analogues"

#### Reference cited in this Technical Information:

- Gjertsen, B.T.; Mellgren, G.; Otten, A.; Maronde, E.; Genieser, H.-G.; Jastorff, B.; Vintermyr, O.K.; McKnight, G.S.; Døskeland, S.O., *J. Biol. Chem.*, **270**, 20599 - 20607 (1995): "Novel (Rp)-cAMPS Analogs as Tools for Inhibition of cAMP-kinase in Cell Culture"