


Technical Information about Rp-8-CPT-cAMPS

Potent membrane permeant, PDE resistant inhibitor of cAMP-dependent protein kinase type I and type II

Update: June 09, 2017 HU



Abbreviation: Rp-8-CPT-cAMPS

Formula	CAS No.	Molecular Weight	UV	BIOLOG Cat. No.
C ₁₆ H ₁₄ ClN ₅ O ₅ PS ₂ ·Na	[129735-01-9]	509.8	λ _{max} 282 nm / ε 16000 / pH7	C 011

Name: 8- (4- Chlorophenylthio)adenosine- 3', 5'- cyclic monophosphorothioate, Rp- isomer (Rp-8-CPT-cAMPS)

Description: Rp-8-CPT-cAMPS is an analogue of the parent second messenger cyclic AMP where the hydrogen in position 8 of the nucleobase is replaced by the lipophilic chlorophenylthio moiety. In addition, the equatorial one of the two exocyclic oxygen atoms in the cyclic phosphate is modified by sulfur (R-isomer). The suffix "p" indicates that R/S nomenclature refers to phosphorus. The compound can be considered as a combination of the well known derivative 8-CPT-cAMP with Rp-cAMPS.

Properties:

- site selective inhibitor of protein kinase A type I and II, preferring site A of type I and site B of type II (Dostmann et al. 1990, Gjertsen et al. 1995),
- metabolic stability towards mammalian cyclic nucleotide-responsive phosphodiesterases,
- high lipophilicity and good membrane permeability while still soluble in aqueous solvents (Genieser 1995),
- can discriminate between protein kinase A (antagonist) and some other cyclic AMP receptors, e.g. ion channels (agonist).

Summing up Rp-8-CPT-cAMPS is a membrane-permeant inhibitor of protein kinase A type I and type II. Its metabolic stability avoids potential side effects through active metabolites. In contrast to common ATP-site inhibitors, Rp-8-CPT-cAMPS occupies cAMP binding sites at the regulatory subunit of PKA and prevents the kinase holoenzyme from dissociation and thus from activation. Due to this working principle, preincubation of the inhibitor prior to the activation step is necessary for optimal results.

Specification: Crystallized or lyophilized sodium salt. Other salt forms of Rp-8-CPT-cAMPS are available upon request. Please keep in mind that equal amounts of the compound may look different in volume depending on humidity. Micromolar quantities are determined by UV at λ_{max}. The agonistic Sp-isomer (Sp-8-CPT-cAMPS) is offered by BIOLOG as well (Cat. No. C 012).


Purity: Typical analysis is better than 99% (HPLC / UV / 282 nm). The product is not sterile and has not been tested for endotoxins.

Stability and Storage: Rp-8-CPT-cAMPS has sufficient stability at room temperature and does not need special care during handling or shipment. Nevertheless, the compound and its solutions should be protected from bright light, stored in the freezer, for longer storage periods preferably in freeze-dried form, since the agonistic 8-CPT-cAMP can be formed slowly by oxidation processes.

Application: Experience shows that applicable concentrations of Rp-8-CPT-cAMPS depend on the type of biosystem, its membrane properties and kinase content. A main application for Rp-8-CPT-cAMPS is to eliminate the first messenger-stimulated phosphorylation by cyclic AMP-dependent protein kinase. For this purpose preincubation (e.g. 20 min) is important since the production of intracellular cyclic AMP initiated by a first messenger is much faster than the antagonist can penetrate the membrane when given extracellularly. Since Rp-8-CPT-cAMPS is hydrolytically stable in mammalian and many other systems there is no danger of degradation during incubation periods. If you have good or moderate results with 8-CPT-cAMP or Rp-cAMPS, you can be sure that Rp-8-CPT-cAMPS will be membrane-permeant in your system as well.

Solubility: Detailed information on the solubility of Rp-8-CPT-cAMPS in water and various buffers are listed in the solubility chart below. Concentrations have been determined at ambient temperature and can be considered as minimum concentrations usually obtainable, however, slight batch-to-batch variations cannot be ruled out. When opening the tube please make sure that no substance is lost within the cap. Please rinse tube walls carefully and preferably use ultrasonic or vortex to achieve total and uniform mixing.

No.	Solvent	Solubility [mM]
I	H ₂ O	100
II	DMSO	100
III	DMF	100
IV	Ethanol 96%	100
V	Methanol	100
VI	PBS, pH 7.4	100
VII	100 mM Na ₂ HPO ₄ , pH 7.0	100
VIII	25 mM Hepes/NaOH, pH 7.2	100
IX	25 mM Tris/HCl, pH 7.4	100



Toxicity and Safety: Since cyclic AMP has multiple tasks in every organism, it is very likely that lipophilic cAMP analogues will interfere with many cell regulation processes *in vivo*. However, due to the rather small quantities to work with no health hazards have been reported. Nevertheless please keep in mind that the *in vivo* properties of this compound are not sufficiently characterized up to now. Avoid skin contact or ingestion and allow only trained personnel to handle the product.

Our products are designed, developed and sold for research purposes only. They are intended for *in vitro* and nonhuman *in vivo* laboratory applications. Any other use requires approval of health authorities.

Not for drug, household or related uses!

Selected References for Rp-8-CPT-cAMPS:

For a detailed list please inquire or visit our website (<http://www.biolog.de>).

Zieba, B.J.; Artamonov, M.V.; Jin, L.; Momotani, K.; Ho, R.; Franke, A.S.; Nepl, R.L.; Stevenson A.S.; Khromov, A.S.; Chrzanowska-Wodnicka, M.; Somlyo, A.V., *J. Biol. Chem.*, **286**, 16681 - 16692 (2011): "The cAMP-Responsive Rap1 Guanine Nucleotide Exchange Factor, EPAC, Induces Smooth Muscle Relaxation by Down Regulation of RhoA Activity"

Roscioni, S.S.; Maarsingh, H.; Elzinga, C.R.S.; Schuur, J.; Menzen, M.; Halayko, A.J.; Meurs, H.; Schmidt, M., *J. Cell. Mol. Med.*, **15**, 1551 - 1563 (2011): "Epac as a Novel Effector of Airway Smooth Muscle Relaxation"

Idevall-Hagren, O.; Barg, S.; Gylfe, E.; Tengholm, A., *J. Biol. Chem.*, **285**, 23007 - 23018 (2010): "cAMP Mediators of Pulsatile Insulin Secretion from Glucose-stimulated Single Beta-Cells"

Grandoch, M.; de Jesus, B.M.; Oude Weernink, P.A.; Weber, A.A.; Jakobs, K.H.; Schmidt, M., *Cell Signal*, **4**, 609 - 621 (2009): "B Cell Receptor-induced Growth Arrest and Apoptosis in WEHI-231 Immature B Lymphoma Cells Involve Cyclic AMP and Epac Proteins"

Waidmann, O.; Pleli, T.; Dvorak, K.; Baehr, C.; Mondorf, U.; Plotz, G.; Biondi, R.M.; Zeuzem, S.; Piiper, A., *J. Biol. Chem.*, **284**, 32256 - 32263 (2009): "Inhibition of the Equilibrative Nucleoside Transporter 1 and Activation of A2A Adenosine Receptors by 8-(4-Chlorophenylthio)-modified cAMP Analogs and their Hydrolytic Products"

Li, S.; Doss, J.C.; Hardee, E.J.; Quock, R.M., *Brain. Res.*, **1038**, 113 - 117 (2005): "Involvement of cyclic GMP-dependent Protein Kinase in Nitrous Oxide-induced Anxiolytic-like Behavior in the Mouse Light/Dark Exploration Test"

Fricke, K.; Schulz, A.; John, H.; Forssmann, W.-G.; Maronde, E., *Endocrinology*, **146**, 2060 - 2068 (2005): "Isolation and Characterization of a Novel POMC-derived Peptide from Hemofiltrate of Chronic Renal Failure Patients"

Rondaj, M.G.; Sellink, E.; Gijzen, K.A.; Ten Klooster, J.P.; Hordijk, P.L.; Van Mourik, J.A.; Voorberg, J., *Arterioscler. Thromb. Vasc. Biol.*, **24**, 1315 - 1320, (2004): "Small GTP-Binding Protein Ral is Involved in cAMP-mediated Release of von Willebrand Factor From Endothelial Cells"

Yang, H.W.; Hu, X.D.; Zhang, H.M.; Xin, W.J.; Li, M.T.; Zhang, T.; Zhou, L.J.; Liu, X.G., *J. Neurophysiol.*, **913**, 1122 - 1133 (2004): "Roles of CaMKII, PKA, and PKC in the Induction and Maintenance of LTP of C-fiber-evoked Field Potentials in Rat Spinal Dorsal Horn"

Waleh, N.; Kajino, H.; Marrache, A.M.; Ginzinger, D.; Roman, C.; Seidner, S.R.; Moss, T.J.M.; Fournon, J.C.; Vazquez-Tello, A.; Chemtob, S.; Clyman, R.I., *Circulation*, **110**, 2326 - 2332 (2004): "Prostaglandin E-2 Mediated Relaxation of the Ductus Arteriosus - Effects of Gestational Age on G Protein-coupled Receptor Expression, Signaling, and Vasomotor Control"

Fricke, K.; Heitland, A.; Maronde, E., *Endocrinology*, **145**, 4940 - 4947 (2004): "Cooperative Activation of Lipolysis by Protein Kinase A and Protein Kinase C Pathways in 3T3-L1 Adipocytes"

Schröder, F.; Klein, G.; Frank, T.; Bastein, M.; Indris, S.; Karck, M.; Drexler, H.; Wollert, K.C., *J. Mol. Cell. Cardiol.*, **37**, 735 - 745 (2004): "Src Family Tyrosine Kinases Inhibit Single L-type: Ca²⁺ Channel Activity in Human Atrial Myocytes"

- Browner, M.C.; Sellak, H.; Lincoln, T.M., *Am. J. Physiol.-Cell Physiol.*, **287**, C88 - C96 (2004): "Downregulation of cGMP-dependent Protein Kinase Expression by Inflammatory Cytokines in Vascular Smooth-muscle Cells"
- Zhu, W.Z.; Wang, S.Q.; Chakir, K.; Yang, D.M.; Zhang, T.; Brown, J.H.; Devic, E.; Kobilka, B.K.; Cheng, H.P.; Xiao, R.P., *J. Clin. Invest.*, **111**, 617 - 625 (2003): "Linkage of β_1 -adrenergic Stimulation to Apoptotic Heart Cell Death through Protein Kinase A-independent Activation of Ca^{2+} /Calmodulin Kinase II"
- Sauzeau, V.; Rolli-Derkinderen, M.; Marionneau, C.; Loirand, G.; Pacaud, P., *J. Biol. Chem.*, **278**, 9472 - 9480 (2003): "RhoA Expression is Controlled by Nitric Oxide through cGMP-dependent Protein Kinase Activation"
- Hsu, K.S.; Huang, C.C.; Liang, Y.C.; Wu, H.M.; Chen, Y.L.; Lo, S.W.; Ho, W.C., *Hippocampus*, **12**, 787 - 802 (2002): "Alterations in the Balance of Protein Kinase and Phosphatase Activities and Age-related Impairments of Synaptic Transmission and Long-Term Potentiation"
- Jang, E.K.; Azzam, J.E.; Dickinson, N.T.; Davidson, M.M.L.; Haslam, R.J., *Brit. J. Haematol.*, **117**, 664 - 675 (2002): "Roles for Both Cyclic GMP and Cyclic AMP in the Inhibition of Collagen-induced Platelet Aggregation by Nitroprusside"
- Christensen, A.E.; de Rooij, J.; Kopperud, R.; Nilsen, O.K.; Dao, K.K.; Schwede, F.; Bos, J.L.; Genieser, H.-G.; Døskeland, S.O., Abstracts of the 13th Protein Kinase Symposium, Oslo/Norway, 2002: "Comparison of the cAMP Binding Sites of EPAC-1 and cAMP-Kinase"
- Gao, Y.S.; Raj, J.U., *Eur. J. Pharmacol.*, **436**, 227 - 233 (2002): "Effects of SQ 22536, an Adenylyl Cyclase Inhibitor, on Isoproterenol-induced Cyclic AMP Elevation and Relaxation in Newborn Ovine Pulmonary Veins"
- Abi-Gerges, N.; Szabo, G.; Otero, A.S.; Fischmeister, R.; Mery, P.F., *J. Physiol.-London*, **540**, 411 - 424 (2002): "NO Donors Potentiate the β -Adrenergic Stimulation of $I_{\text{Ca,L}}$ and the Muscarinic Activation of $I_{\text{K,ACh}}$ in Rat Cardiac Myocytes"
- Han, J.; Kim, N.; Kim, E.; Ho, W.-K.; Earm, Y.E., *J. Biol. Chem.*, **276**, 22140 - 22147 (2001): "Modulation of ATP-sensitive Potassium Channels by cGMP-dependent Protein Kinase in Rabbit Ventricular Myocytes"
- Lubomirov, L.; Gagov, H.; Petkova-Kirova, P.; Duridanova, D.; Kalentchuk, V.U.; Schubert, R., *Br. J. Pharmacol.*, **134**, 1564 - 1570 (2001): "Urocortin Relaxes Rat Tail Arteries by a PKA-mediated Reduction of the Sensitivity of the Contractile Apparatus for Calcium"
- Spicuzza, L.; Belvisi, M.G.; Birrell, M.A.; Barnes, P.J.; Hele, D.J.; Giembycz, M.A., *Br. J. Pharmacol.*, **133**, 1201 - 1212 (2001): "Evidence that the Anti-Spasmogenic Effect of the Beta-Adrenoceptor Agonist, Isoprenaline, on Guinea-Pig Trachealis is not Mediated by Cyclic AMP-dependent Protein Kinase"
- Staples, K.J.; Bergmann, M.; Tomita, K.; Houslay, M.D.; McPhee, I.; Barnes, P.J.; Giembycz, M.A.; Newton, R., *J. Immunol.*, **167**, 2074 - 2080 (2001): "Adenosine 3', 5'-cyclic Monophosphate (cAMP)-dependent Inhibition of IL-5 from Human T Lymphocytes is not Mediated by the cAMP-dependent Protein Kinase A¹"
- Money-Kyrle, A.R.W.; Davies, C.H.; Ranu, H.K.; Ogara, P.; Kent, N.S.; Poole-Wilson, P.A.; Harding, S.E., *Cardiovasc. Res.*, **37**, 532 - 540 (1998): "The Role of cAMP in the Frequency-dependent Changes in Contraction of Guinea-pig Cardiomyocytes"
- Vischer, U.M.; Wollheim, C.B., *Blood*, **91**, 118 - 127 (1998): "Purine Nucleotides Induce Regulated Secretion of von Willebrand Factor - Involvement of Cytosolic Ca^{2+} and Cyclic Adenosine Monophosphate-dependent Signaling in Endothelial Exocytosis"
- Kurtz, A.; Gotz, K.H.; Hamann, M.; Kieninger, M.; Wagner, C., *Amer. J. Physiol.*, **43**, F709 - F717 (1998): "Stimulation of Renin Secretion by NO Donors is Related to the cAMP Pathway"
- O'Neil, R.G.; Leng, L., *Am. J. Physiol.*, **273**, F120 - F128 (1997): "Osmo-Mechanically Sensitive Phosphatidylinositol Signaling Regulates a Ca^{2+} Influx Channel in Renal Epithelial Cells"
- Bonci, A.; Williams, J.T., *J. Neurosci.*, **17**, 796 - 803 (1997): "Increased Probability of GABA Release During Withdrawal from Morphine"
- Maronde, E.; Middendorff, R.; Telgmann, R.; Müller, D.; Hemmings, B.; Tasken, K.; Olcese, J., *J. Neurochem.*, **68**, 770 - 777 (1997): "Melatonin Synthesis in the Bovine Pineal Gland is Regulated by Type II Cyclic AMP-Dependent Protein Kinase"
- Zhou, Y.Y.; Cheng, H.P.; Bogdanov, K.Y.; Hol, C.; Altschuld, R.; Lakatta, E.G.; Xiao, R.P., *Amer. J. Physiol. Heart Circ. Physiol.*, **42**, H1611-H1618 (1997): "Localized cAMP-Dependent Signaling Mediates Beta(2)-Adrenergic Modulation of Cardiac Excitation-Contraction Coupling"
- Schubert, R.; Serebryakov, V.N.; Mewes, H.; Hopp, H.H., *Amer. J. Physiol. - Heart Circ. Phys.*, **41**, H1147-H1156 (1997): "Iloprost Dilates Rat Small Arteries - Role of K-ATP- and K-Ca-Channel Activation by cAMP-Dependent Protein Kinase"
- Stief, C.G.; Ückert, S.; Truss, M.C.; Becker, A.J.; Machtens, S.; Jonas, U., *Urolog. Res.*, **24**, 333 - 337 (1996): "A Possible Role for Nitric Oxide in the Regulation of Human Uretral Smooth Muscle Tone In Vitro"
- Salin, P.A.; Malenka, R.C.; Nicoll, R.A., *Neuron*, **16**, 797 - 803 (1996): "Cyclic AMP Mediates a Presynaptic Form of LTP at Cerebellar Fiber Synapses"
- Azula, F.J.; Alzola, E.S.; Conde, M.; Trueba, M.; Macarulla, J.M.; Marino, A., *Mol. Pharmacol.*, **50**, 367 - 379 (1996): "Thrombin-stimulated Phospholipase C Activity is Inhibited Without Visible Delay by a Rapid Increase in the Cyclic GMP Levels Induced by Sodium Nitroprusside"
- Houillier P.; Chambrey R.; Achard, J.M.; Froissart, M.; Poggioli, J.; Paillard, M., *Kidney Int.*, **50**, 1496 - 1505 (1996): "Signaling Pathways in the Biphasic Effect of Angiotensin II on Apical NA/H Antiport Activity in Proximal Tubule"

Kramer, R.H.; Tibbs, G.R., *J. Neurosci.*, **16**, 1285 - 1293 (1996): "Antagonists of Cyclic Nucleotide-gated Channels and Molecular Mapping of their Site of Action"

Yu, P.Y.; Eisner, G.M.; Yamaguchi, I.; Mouradian, M.M.; Felder, R.A.; Jose, P.A., *J. Biol. Chem.*, **271**, 19503 - 19508 (1996): "Dopamine D-1A Receptor Regulation of Phospholipase C Isoform"

Hemmings, H.C., in *Regulatory Protein Modification: Techniques and Protocols; Neuromethods*, **30**, 121 - 218 (1996); H.C. Hemmings (ed.), Humana Press Inc.: "Protein Kinase and Phosphatase Inhibitors. Applications in Neuroscience"

Boese, M.; Busse, R.; Mulsch, A.; Schini-Kerth, V., *Br. J. Pharmacol.*, **119**, 707 - 715 (1996): "Effect of cyclic GMP-dependent Vasodilators on the Expression of Inducible Nitric Oxide Synthase in Vascular Smooth Muscle Cells: Role of cyclic AMP"

References cited in this Technical Information:

Dostmann, W.R.G.; Taylor, S.S.; Genieser, H.-G.; Jastorff, B.; Døskeland, S.O.; Øgreid, D., *J. Biol. Chem.*, **265**, 10484 - 10491 (1990): "Probing the Cyclic Nucleotide Binding Sites of cAMP-Dependent Protein Kinase I and II With Analogs of Adenosine 3',5'-Cyclic Phosphorothioates"

Gjertsen, B.T.; Mellgren, G.; Otten, A.; Maronde, E.; Genieser, H.-G.; Jastorff, B.; Vintermyr, O.K.; McKnight, G.S.; Døskeland, S.O., *J. Biol. Chem.*, **270**, 20599 - 20607 (1995): "Novel (Rp)-cAMPS Analogs as Tools for Inhibition of cAMP-kinase in Cell Culture"

Genieser, H.-G., Proc. 9. Int. Conf. Second Mess. & Phosphoprot., Nashville, TN, USA, 1995: "Lipophilic Ranking of Biologically Important Cyclic Nucleotides"