


Technical Information about c-diAMP

Update: January 22, 2019 HU



Abbreviation: c-diAMP

Formula	CAS No.	Molecular Weight	UV	BIOLOG Cat. No.
C ₂₀ H ₂₄ N ₁₀ O ₁₂ P ₂ (for free acid)	[54447-84-6]	658.4 (for free acid)	λ _{max} 259 nm / ε 27000 / pH 7	C 088

Name: Cyclic diadenosine monophosphate (c-diAMP) / cyclic bis (3' → 5') diadenylic acid

Description: In c-diAMP two 5'-AMP units are connected to form a cyclic structure.

Properties: c-diAMP was identified as a second messenger that signals DNA integrity in *Bacillus subtilis* during sporulation, and is considered to be an additional purine-based signalling nucleotide with a more general function in prokaryotes (Witte et al. 2008, Oppenheimer-Shaanan et al. 2011). The structurally related bacterial second messenger cyclic diguanosine monophosphate (c-diGMP) is offered as well (Cat. No. C 057).

Specification: Crystallized or lyophilized sodium salt. Please keep in mind that equal amounts of the compound may look different in volume due to sensitivity of the lyophilized form to humidity. The compound can even contract to small volume droplets. Normally the product is located in the conical bottom of the tube. Micromolar quantities are determined by UV at λ_{max}.

Purity: Typical analysis is better than 98% (HPLC / UV / 259 nm). The product is not sterile and has not been tested for endotoxins.

Solubility: c-diAMP is soluble in water (≥ 74 mM, limits have not been determined). Please rinse tube walls carefully and preferably use ultrasonic or vortex to achieve total and uniform mixing. When opening the tube please make sure that no substance is lost within the cap.

Stability and Storage: c-diAMP has sufficient stability at room temperature and does not need special care during handling or shipment. Nevertheless, we recommend that the compound should be stored in the freezer, for longer storage periods preferably in freeze-dried form.

Toxicity and Safety: Please keep in mind, that the *in vivo* properties of this compound are not sufficiently characterized up to now. Avoid skin contact or ingestion and allow only trained personnel to handle the product. Our products are designed, developed and sold for research purposes only! They are intended for *in vitro* and nonhuman *in vivo* laboratory applications. Any other use requires approval of health authorities.

Not for drug, household or related uses!

Selected References for c-diAMP:

Andrade, W.A.; Firon, A.; Schmidt, T.; Hornung, V.; Fitzgerald, K.A.; Kurt-Jones, E.A.; Trieu-Cuot, P.; Golenbock, D.T.; Kaminski, P.-A., *Cell Host Microbe*, **20**, 49-59 (2016): "Group B Streptococcus Degrades Cyclic-di-AMP to Modulate STING-Dependent Type I Interferon Production".

Huynh, T.N.; Luo, S.; Pensinger, D.; Sauer, J.-D.; Tong, L.; Woodward, J.J., *PNAS*, **112**, E747-E756 (2015): "An HD-Domain Phosphodiesterase Mediates Cooperative Hydrolysis of c-di-AMP to Affect Bacterial Growth and Virulence"

Gundlach, J.; Dickmanns, A.; Schröder-Tittmann, K.; Neumann, P.; Kaesler, J.; Kampf, J.; Herzberg, C.; Hammer, E.; Schwede, F.; Kaever, V.; Tittmann, K.; Stülke, J.; Ficner, R., *J.Biol.Chem.*, **290**, 3069 – 3080 (2015): " Identification, Characterization, and Structure Analysis of the Cyclic di-AMP-binding PII-like Signal Transduction Protein DarA"

Sureka, K.; Choi, P.H.; Precit, M.; Delince, M.; Pensinger, D.A.; Huynh, T.N.; Jurado, A.R.; Goo, Y.A.; Sadilek, M.; Iavarone, A.T.; Sauer, J.-D.; Tong, L.; Woodward, J.J., *Cell*, **158**, 1389 - 1401 (2014): "The Cyclic Dinucleotide c-diAMP is an Allosteric Regulator of Metabolic Enzyme Function"

Lolicato, M.; Bucchi, A.; Arrigoni, C.; Zucca, S.; Nardini, M.; Schroeder, I.; Simmons, K.; Aquila, M.; DiFrancesco, D.; Bolognesi, M.; Schwede, F.; Kashin, D.; Fishwick, C.W.; Johnson, A.P.; Thiel, G.; Moroni, A., *Nat. Chem. Biol.*, **10**, 457 - 462 (2014): "Cyclic Dinucleotides bind the C-linker of HCN4 to Control Channel cAMP Responsiveness"

Bai, Y.; Yang, J.; Zarrella, T.M.; Zhang, Y.; Metzger, D.W.; Bai, G., *J. Bacteriol.*, **196**, 614 - 623 (2014): "Cyclic di-AMP Impairs Potassium Uptake Mediated by a c-di-AMP Binding Protein in *Streptococcus pneumoniae*"

Mehne, F.M.P.; Gunka, K.; Eilers, H.; Herzberg, C.; Kaever, V.; Stülke, J., *J. Biol. Chem.* **288**(3),:2004 - 2017 (2013): "Cyclic Di-AMP Homeostasis in *Bacillus subtilis*: Both Lack and High Level Accumulation of the Nucleotide Are Detrimental for Cell Growth"

Yi, G.; Brendel, V.P.; Shu, C.; Li, P.; Palanathan, S.; Kao, C.C., *PLoS One*, **8**, e77846 (2013): "Single Nucleotide Polymorphisms of Human STING Can Affect Innate Immune Response to Cyclic Dinucleotides"

Schaap, P., *IUBMB Life*, **65**, :897 - 903 (2013): "Cyclic Di-Nucleotide Signaling Enters the Eukaryote Domain"

Bai, Y.; Yang, J.; Eisele, L.E.; Underwood, A.J.; Koestler, B.J.; Waters, C.M.; Metzger, D.W.; Bai, G., *J. Bacteriol.*, **195**, 5123 - 5132 (2013): "Two DHH Subfamily Proteins in *Streptococcus pneumoniae* Possess Cyclic Di-AMP Phosphodiesterase Activity and Affect Bacterial Growth and Virulence"

Witte, C.E.; Whiteley, A.T.; Burke, T.P., *mBio*, **4**, e00282-13 (2013): "Cyclic di-AMP is Critical for *Listeria monocytogenes* Growth, Cell Wall Homeostasis, and Establishment of Infection"

Danilchanka O.; Mekalanos J.J., *Cell*, **154**, 962 – 970 (2013): "Cyclic Dinucleotides and the Innate Immune Response"

Corrigan, R.M.; Campeotto, I.; Jeganathan, T.; Roelofs, K.G.; Lee, V.T.; Gründling, A., *Proc. Natl. Acad. Sci. USA*, **110**, 9084 - 9089 (2013): "Systematic Identification of Conserved Bacterial c-di-AMP Receptor Proteins"

Abdul-Sater, A.A.; Grajkowski, A.; Erdjument-Bromage, H.; Plumlee, C.; Levi, A.; Schreiber, M.T.; Lee, C.; Shuman, H.; Beaucage, S.L.; Schindler, C., *Microbes Infect.*, **14**, 188 - 197 (2012): "The Overlapping Host Responses to Bacterial Cyclic Dinucleotides"

Davies, B.W.; Bogard, R.W.; Young, T.S.; Mekalanos, J.J., *Cell*, **149**, 358 - 370 (2012): "Coordinated Regulation of Accessory Genetic Elements Produces Cyclic Di-Nucleotides for *V. Cholerae* Virulence"

Chen, Z.-H.; Schaap, P., *Nature*, **488**, 680 - 683 (2012): "The Prokaryote Messenger c-di-GMP Triggers Stalk Cell Differentiation in *Dictyostelium*"

Corrigan, R.M.; Abbott, J.C.; Burhenne, H.; Kaever, V.; Gründling, A., *PLoS Pathog.*, **7**, e1002217 (2011): "c-di-AMP is a New Second Messenger in *Staphylococcus aureus* with a Role in Controlling Cell Size and Envelope Stress"

Smith, K.D.; Strobel, S.A., *Biochem. Soc. Trans.*, **39**, 647 - 651 (2011): "Interactions of the c-di-GMP Riboswitch with its Second Messenger Ligand"

Gomelsky, M., *Mol. Microbiol.*, **79**, 562 - 565 (2011): "cAMP, c-di-GMP, c-di-AMP and now cGMP: Bacteria use Them All!"

Oppenheimer-Shaanan, Y.; Wexselblatt, E.; Katzhendler, J.; Yavin, E.; Ben-Yehuda, S., *EMBO Rep.*, **12**, 594 - 601 (2011): "c-di-AMP Reports DNA Integrity During Sporulation in *Bacillus subtilis*"

Ebensen, T.; Libanova, R.; Schulze, K.; Yevsa, T.; Morr, M.; Guzmán C.A., *Vaccine*, **29**, 5210 - 5220 (2011): "Bis-(3',5')-cyclic Dimeric Adenosine Monophosphate: Strong Th1/Th2/Th17 Promoting Mucosal Adjuvant"

Abdul-Sater, A.A.; Grajkowski, A.; Erdjument-Bromage, H.; Plumlee, C.; Levi, A.; Schreiber, M.T.; Lee, C.; Shuman, H.; Beaucage, S.L.; Schindler, C., *Microbes Infect.*, **14**, 188 - 197 (2011): "The Overlapping Host Responses to Bacterial Cyclic Dinucleotides"

Tchigvintsev, A.; Xu, X.; Singer, A.; Chang, C.; Brown, G.; Proudfoot, M.; Cui, H.; Flick, R.; Anderson, W.F.; Joachimiak, A.; Galperin, M.Y.; Savchenko, A.; Yakunin, A.F., *J. Mol. Biol.*, **402**, 524 - 538 (2010): "Structural Insight into the Mechanism of c-di-GMP Hydrolysis by EAL Domain Phosphodiesterases"

Woodward, J.J.; Iavarone, A.T.; Portnoy, D.A., *Science*, **328**, 1703 - 1705 (2010): "c-di-AMP Secreted by Intracellular *Listeria monocytogenes* Activates a Host Type I Interferon Response"

Rao, F.; See, R.Y.; Zhang, D.; Toh, D.C.; Ji, Q.; Liang, Z.-X., *J. Biol. Chem.*, **285**, 473 - 482 (2010): "YybT is a Signaling Protein that Contains a Cyclic Dinucleotide Phosphodiesterase Domain and a GGDEF Domain with ATPase Activity"

Pesavento, C.; Hengge, R., *Curr. Opin. Microbiol.*, **12**, 170 - 176 (2009): "Bacterial Nucleotide-based Second Messengers"

Römmling, U., *Sci. Signal.*, **1(33)**, pe39 (2008): "Great Times for Small Molecules: c-di-AMP, a Second Messenger Candidate in Bacteria and Archaea"

Witte, G.; Hartung, S.; Büttner, K.; Hopfner, K.-P., *Mol. Cell*, **30**, 167 - 178 (2008): "Structural Biochemistry of a Bacterial Checkpoint Protein Reveals Diadenylate Cyclase Activity Regulated by DNA Recombination Intermediates"

Simm, R.; Lusch, A.; Kader, A.; Andersson, M.; Römmling, U., *J. Bacteriol.*, **189**, 3613 - 3623 (2007): "Role of EAL-containing Proteins in Multicellular Behavior of *Salmonella enterica* Serovar Typhimurium"